

User Manual for Stainless Steel Smart Watch with 24h HR & Blood Pressure DT92



Specifications

- Display size: 1.28"
- Connection: Bluetooth 4.0
- Waterproof rating: IP67
- Battery capacity: 300mAh
- Working time: 3 - 5 days (depends on user habits)
- Charging time: about 2.5h
- System compatibility: Android, iOS
- Wrist circumference: maximum 200mm, minimum 150mm
- Supported languages: French, English, Portugal, Spain, Italy, Poland, Russian, Arabic, Thai, Chinese, Japanese



1. Charging port
2. Heart rate reading sensor
3. Power On/Off
4. Return Button

WearFit2.0 app installation

Scan the QR code below or search the Google Play Store or App Store to download and install the WearFit2.0 app.



Setting the app for iOS devices: When the device is connected to the watch, the pairing information will pop up. You need to click Agree. App Bluetooth and audio Bluetooth need to be connected at the same time (Bluetooth protocol problem otherwise). The audio Bluetooth can be used only after app Bluetooth is set to connect. Audio Bluetooth will turn off 90 seconds after the watch is turned on, so priority should be given to the connection after the watch is turned on.

Setting the app for Android devices: Enter mobile settings / application management / authorization management / application rights management / find the corresponding app and allow all permissions. If you need to connect the call function and audio Bluetooth at the same time, the watch will turn off audio Bluetooth 90 seconds after powering it on, so it is necessary to pair the watch first thing after powering it on.

WearFit2.0 app introduction

1. On the home page, you can see personal information and additional functions
2. Sleep mode: Record the daily, weekly, and monthly sleep information in the form of a chart.
3. Step mode: You can see the steps, distance traveled and calories burned. It provides you an insight into daily, weekly, and monthly data.
4. Heart rate mode: It provides detailed information on your heart rate on a daily, weekly, and monthly basis. Your heart rate will be read every hour.
5. Blood pressure mode: You get detailed information on your blood pressure on a daily, weekly, and monthly basis. Your blood pressure will be read every hour.
6. Blood oxygen mode: You get detailed information on your blood oxygen on a daily, weekly, and monthly basis. Your blood oxygen will be read every hour.
7. Fatigue mode: Measures your real-time fatigue and displays information hourly.